

HLA-DQA2 antibody



Catalog Number: 111423

Product name

HLA-DQA2 antibody

Specificity

human, mouse; other species not tested.

Antibody description

HLA-DQA2 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, mouse spleen tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 28kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_020056). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:111423(HLA-DQA2 antibody) at dilution of 1:300