

HLA-DRB4 antibody



Catalog Number: 111427

Product name

HLA-DRB4 antibody

Specificity

Human; other species not tested.

Antibody description

HLA-DRB4 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human serum tissue. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of HLA-DRB4 recombinant protein (Accession Number: NM_021983). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111427(HLA-DRB4 antibody) at dilution of 1:600