

# HRASLS5 antibody



Catalog Number: 111459

## Product name

HRASLS5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HRASLS5 Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue, human heart tissue, human ovary tissue, human skin tissue, human spleen tissue. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of HRASLS5 recombinant protein (Accession Number: NM\_001146729). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



Immunohistochemical of paraffin-embedded human placenta using Catalog No:111459(HRASLS5 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:111459(HRASLS5 antibody) at dilution of 1:50 (under 40x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111459(HRASLS5 Antibody) at dilution of

# HRASLS5 antibody



Catalog Number: 111459

---

1:1000