

HNF1A antibody

Catalog Number: 111495

Product name

HNF1A antibody

Specificity

Human, Rat; other species not tested.

Antibody description

HNF1A Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in rat liver tissue, SMMC-7721 cells. Observed molecular weight by Western-blot: 67 kDa

Preparation

This antibody was obtained by immunization of HNF1A recombinant protein (Accession Number: NM_000545). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



rat liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111495(HNF1A antibody) at dilution of 1:1000



IP Result of anti-HNF1A (IP:Catalog No:111495, 4ug; Detection:Catalog No:111495 1:800) with mouse liver tissue lysate 3600ug.