

HNRNPK antibody



Catalog Number: 111509

Product name

HNRNPK antibody

WB: 1:500-1:5000

IP: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

HNRNPK Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Positive FC detected in HepG2 cells. Positive IHC detected in human breast cancer tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 51 kDa

Validations




Preparation

This antibody was obtained by immunization of HNRNPK recombinant protein (Accession Number: XM_017014668). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111509(HNRNPK antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:111509(HNRNPK antibody) at dilution of 1:50 (under 10x lens)

Applications


ELISA, WB, IHC, IF, FC, IP

Dilutions


Recommended Dilution:

HNRNPK antibody

Catalog Number: 111509




Immunofluorescent analysis of HepG2 cells, using HNRNPK antibody Catalog No:111509 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-HNRNPK (IP:Catalog No:111509, 4ug; Detection:Catalog No:111509 1:1000) with HeLa cells lysate 2000ug.



1X10⁶ HepG2 cells were stained with .2ug HNRNPK antibody (Catalog No:111509, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.