

# HORMAD1 antibody

Catalog Number: 111525

## Product name

HORMAD1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HORMAD1 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 45 kDa

## Preparation

This antibody was obtained by immunization of HORMAD1 recombinant protein (Accession Number: XM\_047431808). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:111525(HORMAD1 antibody) at dilution of



1:1000

IP Result of anti-HORMAD1 (IP:Catalog No:111525, 4ug; Detection:Catalog No:111525 1:500) with mouse testis tissue lysate 4000ug.