

# HSD3B7 antibody



Catalog Number: 111553

## Product name

HSD3B7 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

HSD3B7 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue, human kidney tissue, human liver tissue. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 36 kDa

## Preparation

This antibody was obtained by immunization of HSD3B7 recombinant protein (Accession Number: NM\_025193). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111553(HSD3B7 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111553(HSD3B7 antibody) at dilution of 1:50 (under 10x lens)