

HSF4 antibody



Catalog Number: 111560

Product name

HSF4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HSF4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue, mouse brain tissue, mouse lung tissue, rat brain tissue. Observed molecular weight by Western-blot: 62-66 kDa

Preparation

This antibody was obtained by immunization of HSF4 recombinant protein (Accession Number: NM_001538). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111560(HSF4 antibody) at dilution of 1:1000