

ICAM3 antibody



Catalog Number: 111591

Product name

ICAM3 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Immunogen

Specificity

Human; other species not tested.

Antibody description

ICAM3 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, U-937 cells. Positive FC detected in Jurkat cells. Positive IHC detected in human colon tissue. Observed molecular weight by Western-blot: 110-160 kDa, 60 kDa

Preparation

This antibody was obtained by immunization of ICAM3 recombinant protein (Accession Number: NM_002162). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, FC


Dilutions

Recommended Dilution:


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111591(ICAM3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon using Catalog No:111591(ICAM3 antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human colon using Catalog No:111591(ICAM3

ICAM3 antibody



Catalog Number: 111591

antibody) at dilution of 1:100 (under 40x lens)



1X10⁶ Jurkat cells were stained with 0.5ug ICAM3 antibody (Catalog No:111591, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.