

IDI2 antibody

Catalog Number: 111608

Product name

IDI2 antibody

Specificity

human; other species not tested.

Antibody description

IDI2 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 27-30kd

Preparation

This antibody was obtained by immunization of IDI2 recombinant protein (Accession Number: NM_033261). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:111608(IDI2 antibody) at dilution of 1:2000