

IFIH1 antibody

Catalog Number: 111616

General Information

Product name:

IFIH1 antibody

Immunogen:

Applications:

ELISA, WB, IP

Specificity:

Human; other species not tested.

Antibody description:

IFIH1 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 117-125kd

Properties

Preparation:

This antibody was obtained by immunization of IFIH1 recombinant protein (Accession Number: BC111750). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal

Ig Type:Rabbit IgG


Applications

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111616(IFIH1 antibody) at dilution of 1:1500



IP Result of anti-IFIH1 (IP:Catalog No:111616, 4ug; Detection:Catalog No:111616 1:800) with Jurkat cells lysate 2400ug.