

IKIP antibody



Catalog Number: 111646

Product name

IKIP antibody

Immunogen

Specificity

Human, Mouse; other species not tested.

Antibody description

IKIP Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse spleen tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 40kd

Preparation

This antibody was obtained by immunization of IKIP recombinant protein (Accession Number: NM_153687). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:111646(IKIP antibody) at dilution of 1:500