

# IFT80 antibody



Catalog Number: 111671

## Product name

IFT80 antibody

## Immunogen

[Human IFT80 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

IFT80 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, mouse brain tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 88-90 kDa

## Preparation

This antibody was obtained by immunization of IFT80 recombinant protein (Accession Number: NM\_020800). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:111671(IFT80 Antibody) at dilution of 1:1000