

IGBP1 antibody



Catalog Number: 111675

Product name

Polyclonal

IGBP1 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

IGBP1 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, HeLa cells, Jurkat cells, MCF7 cells. Observed molecular weight by Western-blot: 40-45kd

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of IGBP1 recombinant protein (Accession Number: NM_001551). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:111675(IGBP1 antibody) at dilution of 1:400

Clonality