

IgG control antibody

Catalog Number: 111690

Product name

IgG control antibody

Antibody description

IgG control Rabbit Polyclonal antibody. Positive IP detected in Control. Positive WB detected in Protein, Control.

Preparation

This antibody was obtained by immunization of N/A. Purification method: Protein A purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:2000-1:20000

IP: N/A

Validations



IP sample detected with different secondary antibodies. Normal Rabbit IgG (Catalog No:111690) as control.



IP sample detected with different secondary antibodies. Normal Rabbit IgG (Catalog No:111690) as control.