

IGHD antibody

Catalog Number: 111695

Product name

IGHD antibody

Specificity

Human; other species not tested.

Antibody description

IGHD Rabbit Polyclonal antibody. Positive WB detected in Raji cells, Jurkat cells. Positive IHC detected in human tonsillitis tissue. Observed molecular weight by Western-blot: 55-70 kDa

Preparation

This antibody was obtained by immunization of IGHG recombinant protein (Accession Number: BC021276). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:111695(IGHD antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111695(IGHD Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:111695(IGHD Antibody) at dilution of 1:50 (under 40x lens)