

# IGHM antibody

Catalog Number: 111696

---

## Product name

IGHM antibody

## Specificity

Human; other species not tested.

## Antibody description

IGHM Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 75 kDa

## Preparation

This antibody was obtained by immunization of IGHM recombinant protein (Accession Number: BC009851). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:111696(IGHM Antibody) at dilution of 1:300