

# Human IgG antibody



Catalog Number: 111721

## Product name

Human IgG antibody

## Specificity

Human; other species not tested.

## Antibody description

Human IgG Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 25 kDa, 50 kDa, 150 kDa

## Preparation

This antibody was obtained by immunization of Recombinant Protein. Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

## Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:111721(Human-IgG Antibody) at dilution of 1:5000