

# IL36RN antibody



Catalog Number: 111741

## General Information

### Product name:

IL36RN antibody

### Immunogen:

### Applications:

ELISA, WB

### Specificity:

Human; other species not tested.

### Antibody description:

IL36RN Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Properties

### Preparation:

This antibody was obtained by immunization of IL36RN recombinant protein (Accession Number: NM\_012275). Purification method: Antigen affinity purified.

### Formulation:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

### Storage:

Store at -20°C. DO NOT ALIQUOT

### Clonality:

Polyclonal


**Ig Type:**Rabbit IgG

## Applications

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111741(IL1F5 Antibody) at dilution of 1:1000