

# IL36RN antibody



Catalog Number: 111741

## Product name

Polyclonal

IL36RN antibody

## Ig Type

## Specificity

Rabbit IgG

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

IL36RN Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Dilutions

## Preparation

Recommended Dilution:

This antibody was obtained by immunization of IL36RN recombinant protein (Accession Number: NM\_012275). Purification method: Antigen affinity purified.

WB: 1:500-1:5000

## Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Validations

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111741(IL1F5 Antibody) at dilution of 1:1000