

IL3RA antibody



Catalog Number: 111744

Product name

IL3RA antibody

Specificity

Human; other species not tested.

Antibody description

IL3RA Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of IL3RA recombinant protein (Accession Number: NM_002183). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:111744(IL3RA antibody) at dilution of 1:300