

IL-17F antibody

Catalog Number: 111764

Product name

IL-17F antibody

Specificity

Human; other species not tested.

Antibody description

IL-17F Rabbit Polyclonal antibody. Positive WB detected in U-937 cells, K-562 cells. Observed molecular weight by Western-blot: 32-35 kDa

Preparation

This antibody was obtained by immunization of IL-17F recombinant protein (Accession Number: NM_052872). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:111764(IL17F antibody) at dilution of 1:200