

IREB2 antibody

Catalog Number: 111830

Product name

IREB2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

IREB2 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, A549 cells, PC-13 cells. Positive IF detected in A549 cells. Observed molecular weight by Western-blot: 90-105 kDa

Preparation

This antibody was obtained by immunization of IREB2 recombinant protein (Accession Number: NM_001320943). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:111830(IREB2 Antibody) at dilution of 1:300



Immunofluorescent analysis of A549 cells using Catalog No:111830(IREB2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)