

# Integrin beta-1 antibody



Catalog Number: 111862

---

## Product name

Integrin beta-1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Integrin beta-1 Rabbit Polyclonal antibody.  
Positive WB detected in HeLa cells, HT-1080 cells.  
Observed molecular weight by Western-blot: 130 kDa

## Preparation

This antibody was obtained by immunization of Integrin beta-1 recombinant protein (Accession Number: NM\_002211). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111862(ITGB1 antibody) at dilution of 1:300