

INTS3 antibody

Catalog Number: 111872

Product name

INTS3 antibody

Immunogen

[Human INTS3 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

INTS3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, mouse brain tissue, mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 118kd

Preparation

This antibody was obtained by immunization of INTS3 recombinant protein (Accession Number: BC025254). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111872(INTS3 antibody) at dilution of 1:500



IP Result of anti-INTS3 (IP:Catalog No:111872, 4ug; Detection:Catalog No:111872 1:500) with mouse testis tissue lysate 4000ug.