

# JAKMIP2 antibody



Catalog Number: 111881

## Product name

Polyclonal

JAKMIP2 antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

JAKMIP2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 95kd

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of JAKMIP2 recombinant protein (Accession Number: NM\_014790). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111881(JAKMIP2 antibody) at dilution of 1:1500