

KANK2 antibody



Catalog Number: 111909

Product name

KANK2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KANK2 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse lung tissue, mouse heart tissue, mouse kidney tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 95kd

Preparation

This antibody was obtained by immunization of KANK2 recombinant protein (Accession Number: NM_001379563). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:111909(KANK2 antibody) at dilution of 1:800



Immunofluorescent analysis of HepG2 cells, using KANK2 antibody Catalog No:111909 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).