

# KBTBD5 antibody



Catalog Number: 111921

## Product name

KBTBD5 antibody

## Specificity

Human; other species not tested.

## Antibody description

KBTBD5 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 69 kDa

## Preparation

This antibody was obtained by immunization of KBTBD5 recombinant protein (Accession Number: NM\_152393). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:111921(KBTBD5 Antibody) at dilution of 1:300