

KCNB1-Specific antibody



Catalog Number: 111926

Product name

KCNB1-Specific antibody

Specificity

Human; other species not tested.

Antibody description

KCNB1-Specific Rabbit Polyclonal antibody.
Positive IF detected in SH-SY5Y cells, HEK-293 cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004975).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:111926(KCNB1-Specific Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG