

# ITPKA antibody



Catalog Number: 111957

## Product name

ITPKA antibody

Recommended Dilution:

WB: 1:500-1:5000

## Immunogen

[Human ITPKA Recombinant protein \(GST tag\)](#)

IP: 1:200-1:1000

IHC: 1:20-1:200

## Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:20-1:200

## Validations

### Antibody description

ITPKA Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Positive IF detected in HepG2 cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 51 kDa




### Preparation

This antibody was obtained by immunization of ITPKA recombinant protein (Accession Number: NM\_002220). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111957(ITPKA antibody) at dilution of 1:500

### Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:111957(ITPKA antibody) at dilution of 1:100 (under 10x lens)

### Storage

Store at -20°C. DO NOT ALIQUOT

### Clonality

Polyclonal

### Ig Type

Rabbit IgG


### Applications

ELISA, WB, IHC, IF, IP


### Dilutions

# ITPKA antibody


Catalog Number: 111957



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:111957(ITPKA antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using ITPKA antibody Catalog No:111957 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent



DNA dye).

IP Result of anti-ITPKA (IP:Catalog No:111957, 4ug; Detection:Catalog No:111957 1:400) with mouse brain tissue lysate 4000ug.