

# KIAA1217 antibody



Catalog Number: 111991

## Product name

KIAA1217 antibody

## Immunogen

[Human KIAA1217 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

KIAA1217 Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, PC-3 cells. Observed molecular weight by Western-blot: 100-160 kDa

## Preparation

This antibody was obtained by immunization of KIAA1217 recombinant protein (Accession Number: BC018764). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:111991(KIAA1217 Antibody) at dilution of 1:600