

KIAA0087 antibody



Catalog Number: 112042

Product name

KIAA0087 antibody

Specificity

Human; other species not tested.

Antibody description

KIAA0087 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, human brain tissue. Positive WB detected in fetal human brain tissue, SW 1990 cells. Observed molecular weight by Western-blot: ~40 kDa

Preparation

This antibody was obtained by immunization of KIAA0087 recombinant protein (Accession Number: BC128257). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:200-1:2000


IHC: 1:20-1:200

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112042(KIAA0087 Antibody) at dilution of 1:600

Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:112042(KIAA0087 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:112042(KIAA0087 Antibody) at dilution of 1:50 (under 40x lens)