

# KIR2DS4 antibody



Catalog Number: 112061

## Product name

KIR2DS4 antibody

## Specificity

Human; other species not tested.

## Antibody description

KIR2DS4 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, Jurkat cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 33 kDa

## Preparation

This antibody was obtained by immunization of KIR2DS4 recombinant protein (Accession Number: NM\_012314). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:112061(KIR2DS4 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using KIR2DS4 antibody Catalog No:112061 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).