

# LIAS antibody

Catalog Number: 112220

## Product name

LIAS antibody

## Immunogen

[Human LIAS Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

LIAS Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 42kd

## Preparation

This antibody was obtained by immunization of LIAS recombinant protein (Accession Number: NM\_006859). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112220(LIAS antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:112220(LIAS antibody) at dilution of 1:100 (under 40x lens)