

# LIPI antibody



Catalog Number: 112244

---

## Product name

LIPI antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

LIPI Rabbit Polyclonal antibody. Positive WB detected in A549 cells, mouse lung tissue. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of LIPI recombinant protein (Accession Number: BC148691). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:112244(LIPI Antibody) at dilution of 1:300