

# Kv4.2 antibody



Catalog Number: 112259

## General Information

### Product name:

Kv4.2 antibody

### Immunogen:

### Applications:

ELISA, WB

### Specificity:

Human, Mouse, Rat; other species not tested.

### Antibody description:

Kv4.2 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, DU 145 cells, HeLa cells. Observed molecular weight by Western-blot: 70-80 kDa

## Properties

### Preparation:

This antibody was obtained by immunization of Kv4.2 recombinant protein (Accession Number: BC110449). Purification method: Antigen affinity purified.

### Formulation:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

### Storage:

Store at -20°C. DO NOT ALIQUOT

### Clonality:

Polyclonal


**Ig Type:** Rabbit IgG

## Applications

Recommended Dilution:

WB: 1:500-1:5000

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:112259(KCND2 antibody) at dilution of 1:300