

# LAMP1 antibody

Catalog Number: 112275

---

## Product name

LAMP1 antibody

## Specificity

Human; other species not tested.

## Antibody description

LAMP1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, human brain tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 110 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_005561). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112275(LAMP1 Antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells using Catalog No:112275(LAMP1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)