

LNPEP antibody

Catalog Number: 112289

Product name

LNPEP antibody

Specificity

Human; other species not tested.

Antibody description

LNPEP Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, COLO 320 cells, human heart tissue. Positive IHC detected in human heart tissue, human kidney tissue. Observed molecular weight by Western-blot: 117-125 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005575). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112289(LNPEP antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human heart using Catalog No:112289(LNPEP antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:112289(LNPEP antibody) at dilution of 1:50 (under 40x lens)