

# L-Plastin antibody



Catalog Number: 112307

## Product name

L-Plastin antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

L-Plastin Rabbit Polyclonal antibody. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 64-70 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_002298). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:112307(LCP1 antibody) at dilution of 1:3000