

# LPPR2 antibody



Catalog Number: 112310

---

## Product name

LPPR2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

LPPR2 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle tissue. Positive WB detected in 3T3-L1 cells, HeLa cells. Observed molecular weight by Western-blot: 37 kDa

## Preparation

This antibody was obtained by immunization of LPPR2 recombinant protein (Accession Number: NM\_022737). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:50-1:500

## Validations



3T3-L1 cells were subjected to SDS PAGE followed by western blot with Catalog No:112310(LPPR2 Antibody) at dilution of 1:1000