

LSP1 antibody



Catalog Number: 112353

Product name

LSP1 antibody

Specificity

human; other species not tested.

Antibody description

LSP1 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, U-937 cells. Observed molecular weight by Western-blot: 52kd

Preparation

This antibody was obtained by immunization of LSP1 recombinant protein (Accession Number: NM_002339). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:112353(LSP1 antibody) at dilution of 1:500