

LTB4R antibody

Catalog Number: 112358

Product name

LTB4R antibody

Specificity

Human; other species not tested.

Antibody description

LTB4R Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, human breast cancer tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_181657).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:112358(LTB4R Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:112358(LTB4R Antibody) at dilution of 1:50 (under 40x lens)