

MAEA antibody



Catalog Number: 112371

Product name

MAEA antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MAEA Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse pancreas tissue. Positive IHC detected in human breast cancer. Observed molecular weight by Western-blot: 33kd,36kd

Preparation

This antibody was obtained by immunization of MAEA recombinant protein (Accession Number: NM_001017405). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112371(MAEA antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:112371(MAEA antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:112371(MAEA antibody) at dilution of 1:50

MAEA antibody



Catalog Number: 112371

(under 40x lens)