

MAGED4 antibody



Catalog Number: 112392

Product name

MAGED4 antibody

Specificity

human; other species not tested.

Antibody description

MAGED4 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, SH-SY5Y cells. Observed molecular weight by Western-blot: 42-60kd

Preparation

This antibody was obtained by immunization of MAGED4 recombinant protein (Accession Number: NM_030801). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:112392(MAGED4 antibody) at dilution of 1:1000