

LY6G5B antibody

Catalog Number: 112400

Product name

LY6G5B antibody

Specificity

Human; other species not tested.

Antibody description

LY6G5B Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human colon cancer tissue.

Preparation

This antibody was obtained by immunization of LY6G5B recombinant protein (Accession Number: NM_021221). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:112400(LY6G5B Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:112400(LY6G5B Antibody) at dilution of 1:50 (under 40x lens)