

Marcks antibody

Catalog Number: 112449

Product name

Marcks antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Marcks Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, mouse pancreas tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 88kd

Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM_008538). Purification method: Protein A purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112449(Marcks antibody) at dilution of 1:2000



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:112449(Marcks antibody) at dilution of 1:200 (under 10x lens)