

MAGI1 antibody



Catalog Number: 112454

Product name

MAGI1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

MAGI1 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, human brain tissue. Observed molecular weight by Western-blot: 110-120kd,165-170kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004742). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:112454(MAGI1 antibody) at dilution of 1:500