

MAVS; VISA antibody



Catalog Number: 112499

Product name

MAVS; VISA antibody

Immunogen

Specificity

Human; other species not tested.

Antibody description

MAVS; VISA Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, HeLa cells, human heart tissue, Jurkat cells. Positive IF detected in HeLa cells. Positive IHC detected in human skin tissue. Observed molecular weight by Western-blot: 70-75kd

Preparation

This antibody was obtained by immunization of MAVS; VISA recombinant protein (Accession Number: NM_020746). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IP, IHC

Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112499(MAVS; VISA antibody) at dilution of



Immunofluorescent analysis of HeLa cells, using MAVS antibody Catalog No:112499 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

MAVS; VISA antibody


Catalog Number: 112499



IgG (red).

IP Result of anti-MAVS; VISA (IP:Catalog No:112499, 3ug; Detection:Catalog No:112499 1:1000) with HEK-293 cells lysate 1700ug.

Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:112499(MAVS; VISA Antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:112499(MAVS; VISA Antibody) at dilution of 1:200 (under 40x lens)

