

# MCFD2 antibody



Catalog Number: 112546

---

## Product name

MCFD2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MCFD2 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells, HeLa cells.

## Preparation

This antibody was obtained by immunization of MCFD2 recombinant protein (Accession Number: NM\_139279). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IF

## Dilutions

Recommended Dilution:

IF: 1:20-1:200

## Validations



Immunofluorescent analysis of HepG2 cells using Catalog No:112546(MCFD2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)