

MICA antibody



Catalog Number: 112608

Product name

MICA antibody

Specificity

human; other species not tested.

Antibody description

MICA Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells, SH-SY5Y cells. Observed molecular weight by Western-blot: 43kd,50-55kd

Preparation

This antibody was obtained by immunization of MICA recombinant protein (Accession Number: BC016929). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112608(MICA antibody) at dilution of 1:500